

CELL SMEAR CENTRIFUGE

Features

- Drive directly by large torque DC motor, with the speed range of 100~4000rpm, microprocessor control and simple operation.
- The protection from over-speed, imbalance and door-cover, which can ensure the safe operation of the instrument.
- It is consistent with CE security certification, and passed ISO9001 and ISO13485 quality system certification.



Technical Parameters

Model	CCL-4ZXB
Electric Supply	AC220V & 110V/50-60Hz
Speed Accuracy	±30rpm
Power	90W
Max Capacity	12X0.5ml
Max Speed	3000r/min
Max RCF	2300xg
Interference Suppression	En61326
Dimension	275x390x300mm(LxWxH)
Weight	20kg
Noise level at max speed	≤ 50dB (A)

Rotor

Number	Max Speed	Max RCF	Capacity
No.1angle Rotor	3000r/min	730xg	6x0.5ml
No.2 Angle Rotor	3000r/min	2300xg	12x0.5ml