



- Easy programming to achieve high accuracy of speed control
- Easy operation with intuitive menu and easy-to-read display
- Rotor speed can be set and displayed by RPM or G-force
- Maintenance-free brushless DC motors
- Quick spin function with press and hold pulse key
- Lid with electric lock and automatic release function saves processing time and prevents samples overheating
- Automatic self-diagnosis on startup showing accumulative running time and last running parameters
- Braking time 2 gears adjustable(for A6-50P rotor)
- A12-10P and A6-50P are not interchangeable

Specifications

CCL-DM0412

Speed Range	300-4500rpm,increment: 100rpm
Max. RCF	2490×g
Speed accuracy	±20rpm
Rotor Capacity	1 type
Run Time	30sec-99min/Continuous
Motor	Brushless DC moto
Display	LCD
Safety Devices	Door interlock,Overspeed detector; Automatic internal diagnosis
Acceleration /Braking time [Sec]	20s↑13s↓
Power	Single-phase, 110V-240V,50Hz/60Hz, 3A 70W
Dimension [D×W×H]	301×354×217mm
Weight	6kg
Certification	CE, UL, IVD
Advanced Features	Speed/RCF switch; Short-time run function; sound-alert function

China Care Medical Equipment Co.,Ltd
 Adress:Room B101, Zhongbang Business Building, PuXing
 Road No. 87, Tianhe, Guangzhou,
 Chinaweb:www.chinacaremedical.com