

High Speed Micro-Centrifuges

D3024 High Speed Micro Centrifuge



2-point lid-lock mechanism – ensures secure closure
Automatic lid opening when run is complete

High strength aluminum alloy rotor can be fully autoclaved with unlimited service life

Special multi-flow air cooling design (Patent No. 201120322916.0) helps keep lower rotor temperature

CE FC

Intuitive symbols – display run status

Full-digit speed display (rpm/Rcf) – easy identification of parameters

Centrifuge goes into a self diagnostics when switched on and then displays accumulative running time and last running parameters

The short spin function is used for quick spins by pressing and holding PULSE key



Special multi-flow air cooling design protects sensitive samples (D3024)
FastTemp function for fast pre-cooling and maximum temperature accuracy inside the rotor (D3024R)



China Care Medical Equipment Co.,Ltd
Address:Room B101,
Zhongbang Business Building,
PuXing Road No. 87, Tianhe,
Guangzhou,
Chinaweb:www.
chinacaremedical.com

