


Marie Stopes® MVA Syringe and Cannula Product Brief

Safe · discreet · effective · client-friendly



China Care Medical Equipment Co.,Ltd
Address:Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
China web:www.chinacaremedical.com





MANUAL VACUUM ASPIRATION (MVA)

- A safe and effective method of surgical abortions for up to 12 weeks of pregnancy

Discreet and Client-Friendly

- Quiet procedure meets clients' needs for privacy
- Enables faster recovery times compared to operating room procedures


Safe and Effective

- Can be performed without the need for a specialist set up or operating theatre
- 98% effective with low complication rates, recommended by the WHO
- Safe alternative to sharp curettage

Affordable

- Reduces procedural costs by performing procedures in an outpatient setting
- No need for general anesthesia

Single use only
CFDA approved
CE-marked
and ISO 13485 compliant



Marie Stopes® Single-valve Syringe

24.5cm length, 3cm diameter, 60ml volume, 375 mm Hg vacuum

Easy to operate

Marie Stopes® MVA Flexible Cannula

Cannula is named by its diameter (e.g. 4mm cannula's diameter is 4mm).

Current cannula sizes are 4mm, 5mm, 6mm, 7mm, 8mm, 9mm, 10mm, 11mm and 12mm.

23.0cm length


Marie Stopes International (MSI)



Marie Stopes International (MSI) is a leading provider of family planning and reproductive healthcare globally. It works in 37 countries around the world delivering services to over 20 million women and men every year to ensure the rights of women to have children by choice, not chance.

As part of reproductive health services, MSI provides abortions using the safe and simple technique which called manual vacuum aspiration (MVA).

MS HEALTH PTY LTD is MSI's wholly owned subsidiary in China and is the exclusive distribution agent in China and overseas for Marie Stopes® MVA product line.



Product Portfolio

- Single-valve Syringe
- MVA Flexible Cannula
(Size 4mm/5mm/6mm/7mm/8mm/9mm/10mm/11mm and 12mm)
- Cervical Probe and Dilator
- MVA Kits
(1 syringe and 2 cannulas)