



Electro-thermo dry box

Electrothermal thermostatic incubator is a necessary equipment for fungi storage and biology culture in scientific research applications. It applies microprocessor, intelligent temperature controller, and durable glass door design, temperature control and easy for observation.

model	power supply	temperature range	temp accuracy	power	work room
CCL-P25	220V50HZ	RT+5℃-65℃	±1℃	150W	30 X 30 X 28
CCL-P50	220V50HZ	RT+5℃-65℃	±1℃	250W	36 X 36 X 42
CCL-P80	220V50HZ	RT+5℃-65℃	±1℃	300W	42 X 42 X 50
CCL-P160	220V50HZ	RT+5℃-65℃	±1℃	400W	50 X 50 X 65
CCL-P270	220V50HZ	RT+5℃-65℃	±1℃	600W	60 X 60 X 75