

Characters


- The software operation is simple and audiovisual with Chinese manifestation.
- The software has strongly analytical sketch function with a high-speed printer.
- Operators compile program arbitrarily.
- The bar code scanner can read directly sample ID number.

Specification

Optics	Top-class flat panel scanner	Scan speed	24 specimen/60seconds
Scan length	216×350mm	Scan way	Automatic, handcraft, constant regulation and enlargement
Operation system	Windows XP	Area name	Several areas can be set
Power	220V 50Hz		
Edit characteristics	The summit of the configure is automatic and handcraft analysis. The baseline can be adjusted. The district can be deleted. The various report format can be set.		
Disease case data	10 compiled contents, such as disease member name, age, sex, time, department etc.		
Ratio value	The rate of the different or some areas can be established.		
Connect	USB RS-232, ASC II and STM output connect		







Data processing system



Operation interface



EPSON Scanner



China Care Medical Equipment Co., Ltd

www.chinacaremedical.com


Cellphone/Whatsapp: 008613926480041 Email: info@chinacaremedical.com

Auto Agarose Gel Electrophoresis System



Features

- ★Auto sample, application, electrophoresis
- ★The time of sample, application, electrophoresis can be set
- ★Sample Blade applicator is clear and high solution
- ★Adjustable temperature of electrophoresis slot
- ★Reagent specification: 4, 6, 12, 24, 30, 36 test/kits
- ★Output voltage: 60-500V
- ★Output current: 0-100mA



Parameters

- ★Hb
- ★SPE
- ★UPE
- ★HbA1c
- ★HDL, VLDL, LDL
- ★ALK-P
- ★Lipo
- ★CK

