

CCE-3312B 12ch Product Introduction

Product Photos :



Main Unit



Accessories



Features :

- Portable delicate design, big color LCD to display and soft alphanumeric keyboard operation
- Accurate pulse pace identification and excellent auto-interpretation function
- High accuracy digital filter, auto baseline adjustment
- 12 leads synchronously acquisition and display; Recording mode : 3ch+++, 6ch, 6ch+, 12ch
- Working modes: Manual / Auto / Analysis / Storage
- 800×480 graphic 7 inch color LCD to display detailed ECG information
- 250 patient cases save and replay (SD card storage is optional)
- Built-in rechargeable Li-ion battery supports continuous work for about 4 hours.
- 215mm*20m roll or 210*140mm folder, high resolution thermal paper

- Rs232 and USB interface (USB disk storage,laser printer printing , PC ECG software are optional)
- Language : English, Chinese, French, Italian, Spanish, Russian
- 300mm*230mm*75mm, 2.8kgs

Technical Specification :

Lead	Standard 12 leads	Time constant	$\geq 3.2s$
Lead acquisition	12bit/1000Hz (synchronously 12 leads)	Frequency response	0.05~160Hz(-3db)
Work mode	Man / Auto/Analysis/Storage	Noise level	$< 15 \mu V_{p-p}$
Filter	AC Filter : 50Hz / 60Hz EMG Filter : 25Hz / 45Hz Anti-Drift Filter : 0.15Hz (Adaptive)	Inter-channel Interference	$\leq 0.5mm$
CMMR	$\geq 120dB$ (with AC FILTER)	Sensitivity	2.5, 5, 10, 20, 40, Auto (mm/mV)
Input circuit	Floating; Protection circuit against Defibrillator effect	Recording mode	3ch+++, 6ch, 6ch+, 12ch Default Man. Mode is 3ch+++ format Default Auto. Mode is 12ch format
Input Impedance	$\geq 50M\Omega$	Paper speed	6.25、12.5、25、50mm/s ($\pm 3\%$)
Input circuit current	$\leq 0.05 \mu A$	Paper	215mm *20m Roll or 210*140mm folder paper
Patient current leakage	$< 10 \mu A$	Communication	RS232/ USB port (optional)
Calibrating voltage	$1mV \pm 2\%$	Power Supply	AC 110-230V($\pm 10\%$),50/60Hz ($\pm 1Hz$),60VA DC Built-in rechargeable battery,14.4V(4400mAh)
Voltage tolerance	$\geq \pm 500mV$	Dimension/Weight	347mm $\times$ 293mm $\times$ 83mm/4.8Kg

Package :

1 unit / small carton : 43cm*37cm*22cm, net weight 4.8kgs, gross weight 7kgs,

4 units / big carton : 75cm*43cm*46cm, net weight 20kgs, gross weight 28kgs,



China Care Medical Equipment Co.,Ltd
Address:Room B101, Zhongbang Business
Building, PuXing Road No. 87, Tianhe,
Guangzhou,
China web: www.chinacaremedical.com