



CCS-BS50

Max Output Power: 30W

1 Characteristic

- Control the power by using the power adjustment button of the ESU pencil so as to control output power and the operation time is shortened.
- Dual microprocessors (one is to control function & the second is to terminate activation if a hazard or safety issue is detected).
- Memory Function: Set up the requested function parameter after turn-off. At the next startup the previous setting parameter will be displayed.
- From 10W to 30W, the output of each watts can be manipulated precisely. And below 10W, the adjustment is done at the 0.1W amplitude by up and down, all of micro surgery is effected by this machine.
- Stimulate the argon system for activating the spray coagulation.
- Can be used with disposable electrosurgical pad.



2 Output Mode

Mono-polar Coagulation: (30W, Load: 500 Ω) Bipolar Coagulation 1: (30W, Load: 100 Ω)

3 Scope of Application

Dermatology, Plastic Surgery, Proctology, Gynecology, Ophthalmology, Urology, Otolaryngology, Dental/Oral Surger

4 Technical Parameter:

Operating ambient temperature	5°C~40°C
Relative humidity range	$\leq 80\%$
Atmosphere pressure range	86.0 ~ 106.0kPa
Power Supply	220/110V, 50/60Hz
Operating Frequency	450kHz $\pm$ 50kHz

5 Transportation Data

Size: 22.5cm (L) *19.5cm (W) *12cm (H) (ESU Generator size) ;
33 (L) *25 (W) *26.5cm (H) (outer carton)
Weight: 2.5kg (Net Weight) ; 4 kg (Gross Weight)

6 Standard Matching Accessories

Disposable power control pencil (optional):2pcs
Or reusable power control pencil (optional):1pc
Disposable solid electrosurgical pad and reusable cable (optional): 1pc
Or reusable stainless steel pad and cable (optional)
Power Cord: 1pc
Reusable Bipolar Forceps and Cable: 1pc
Individual Pedal Switch: 1pc
Instruction Manual: 1pc

China Care Medical Equipment Co.,Ltd
Adress:Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
Chinaweb:www.chinacaremedical.com