

Vertical System



Features

- Easy to cast the gel
- Transparent, no obstructions to visibility
- Distinctive design for elimination of heat, there is a room between the gel room and the main tank plate, which is for buffer solution and buffer solution will eliminate the heat produced during the running.
- Simple and easy installation of the tank.
- High quality and clear electrophoresis patterns can be obtained

Specifications

- Gel dimensions: 480 X 170 (mm)
- Runs up to 26 or 32 samples
- Comb thickness: 0.4 (mm)
- Buffer volume: about 750 ml
- Size (L X W X D): 235 X 180 X 550 (mm); weight: 9 kg
- Recommended Power Supplies: DYY-10C, DYY-12C

Applications

For DNA Sequencing, Differential display, AFLP or SSCP research and Oligonucleotide analysis etc.



China Care Medical Equipment Co.,Ltd
Address: Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
China web: www.chinacaremedical.com