


01 CCL-MNES Mini-Electrophoresis System



Specification

Input Power	AC100~120V, 50~60Hz / AC200~240V, 50~60Hz
Output Power	DC35V / DC50V / DC100V
Bath Dimensions (mm)	110 x 121 x 43mm
Volume of Bath	230ml
Construction of Bath	PC+ABS with high temperature resistance
Timer range	1min~99min
Max.Power	40W
Dimensions (mm)	190 x 130 x 60mm

Accessories

Tooth wide (width x thickness x no.)	3mm x 1mm x 22 / 5.6mm x 1mm x 12 / 3mm x 1mm x 18 / 5.6mm x 1mm x 10
Gel Dimension (W x D)mm	126 x 126mm-1pcs
Gel Migration Board (W x D)mm	110 x 60mm-2pcs / 54 x 60mm-4pcs

Ordering information

No.	Code	Description
01	AS-14010-00	Mini-ES1, AC100~120V, 50~60Hz
02	AS-14020-00	Mini-ES2, AC200~240V, 50~60Hz

02 K37-24 Gel Card Incubator

Description

The K37-24 Gel card incubator can be suitable for many types Gel cards from different manufacturers

Features

1. 2 independent incubation zones with a capacity for 12 Gel cards
2. Incubator for 24-Gel Cards in racks at 37°C
3. Digital interface screen for control and visualization of remaining incubation time and real temperature for each zone



China Care Medical Equipment Co.,Ltd
 Address:Room B101, Zhongbang Business Building, PuXing
 Road No. 87, Tianhe, Guangzhou,
 Chinaweb:www.chinacaremedical.com