


MEDICAL

Anesthesia consumable for endotracheal intubation



Perfect ergonomic design



Unique dermatoglyph anti-slip design,
increase the friction in hand



Highlight LED cold light source,
the light is better, tube never get hot




High button & lithium battery,
use more secure

One-piece Design Disposable Laryngoscope



ISO13485

ISO9001



- The UShare one piece blade and handle combination is a versatile, high quality, airway management tool.
- This disposable instrument offers the features necessary to assist in quality intubation.



Product Description

Model	Size	Blade Length	Handle cover color
CCEL01	Large	159.7mm	Purple
CCEL02	Medium	139.6mm	Blue
CCEL03	Small	116.1mm	Green
CCEL04	Pediatric	107.8mm	Pink
CCEL05	Neonatal	90mm	Yellow

Product Features

Advise: the warnings and notice are listed in User Instruction.

- Unique unibody design, maximum guarantee stable operation, save operating time.
- Super high performance medical ABS material, large bearing capacity.
- Finely finished, no-traumatic edges with a real sooth tip.
- Independent sterile packaging and disposable use prevent cross infection.
- Five different models to meet different patients needs .
- No cleaning and disinfection required, discarded after use, the use is more convenient.
- Endurable LED light source, up to 3 hours of constant light quality.



China Care Medical Equipment Co.,Ltd
 Address: Room B101, Zhongbang Business Building,
 PuXing Road No. 87, Tianhe, Guangzhou,
 Chinaweb: www.chinacaremedical.com