



Video Gastroscope

CMOS 1,000,000 pixels HD resolution



Technical Specification

Item	Gastroscope
Diameter of Distal-end	Ø9.2mm/9.8mm
Diameter of Insert Tube	Ø9.0mm/9.6mm
Biopsy channel	Ø 2.8mm
Working Length	1030mm
Total Length	1330mm
View of Field	140°
Depth of View	3-100mm
Resolution	1,000,000 CMOS
Tip Deflection	U:210° D:90° L/R:100°
Remark: We can provide OEM service , technical details can be customized.	

China Care Medical Equipment Co.,Ltd
Adress:Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
Chinaweb:www.chinacaremedical.com