



Handheld Digital Medical Scope Series

Product Features:

- LCD screen: 3.5 inch
- Screen can be adjusted to 360° (left to right) ,120° (front to back)
- 24 pcs disposable Len covers in standard configuration
- 3.7 V 18650 lithium battery can be removable.

Product specifications:

- Resolution: 720 x 480
- Light source: Natural white (LED)
- LCD screen: 3.5 inch high-definition LCD
- Image format: JPEG, BMP
- Media: Micro SD flash card
- Power supply: rechargeable lithium battery
- External power: 110~220 VAC, 50/60 Hz

Product Appearance:



China Care Medical Equipment Co.,Ltd
 Address:Room B101, Zhongbang Business Building, PuXing
 Road No. 87, Tianhe, Guangzhou,
 Chinaweb:www.chinacaremedical.com