



Handheld Digital Medical Scope Series



Quality, Professional, Innovational

China Care Medical Equipment Co.,Ltd
Adress:Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
Chinaweb:www.chinacaremedical.com

HANDHELD DIGITAL MEDICAL SCOPE SERIES

Product Features

- Specially designed for the department of ENT. It can be used to examine eardrum, nasal cavity and throat with different camera probe.
- Avoid high cost of the traditional ENT examination equipment, such as light source, CCD, Monitor etc. All the functions are in one handheld device.
- Light weight and portable, suitable for outpatient department, emergency room and medical ward.
- High resolution image. The SD card can store multiple frozen images and videos, and can be connected to computer to print patient reports.
- 3 inches high definition LCD screen, watching and taking photos at the same time.
- It can be connected to larger size monitor, synchronized to the external bigger screen.

Product Appearance




Original Images



normal ear drum



normal ear drum



ear drum with otitis media

Product Specifications

- Charge: 110V~220V AC 50/60HZ Output DC 5V 500MA
- LCD: 3 inches special LCD screen, minimum color difference for different viewing angle.

Applicable Camera Probes

- The device is applicable for viewable examination for ear, nose and throat with different camera.



otoscope(standard)



rhinoscope(optional)



laryngoscope(optional)

China Care Medical Equipment Co.,Ltd
Adress:Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
Chinaweb:www.chinacaremedical.com