

CCEU32 Reusable Uretero-
Renoscope **Technical Specifications**

Optical System	Pixels: 160K
	Field of view: 120°
	Depth of view: 3-100mm
	Illumination: optical fiber
Insertion Tube	Maximum outer diameter: 8.5fr
	Working length: 650mm
Instrument Channel	Working channel: 3.6fr
Bending Section	Angulation range: 285°
Function	Photo
	Video
	White Balance
	Storage: 16G
Connections	DVI
	S-Video
	CVBS



China Care Medical Equipment Co.,Ltd
Address:Room B101, Zhongbang Business Building,
PuXing Road No. 87, Tianhe, Guangzhou,
China web:www.chinacaremedical.com

**AN ENDOSCOPE
SOLUTION
COMPANY**

**Reusable
Uretero-Renoscope
Always Available
Cost Effective**

China Care Medical Equipment Co.,Ltd

Cellphone/Whatsapp: 008613926480041 Skype: chinacaremed Wechat:
Joseph20130708

Facebook/Youtube/Instragm/Linkedln: china care medical

Website: www.chinacaremedical.com Email: info@chinacaremedical.com

www.chinacaremedical.com






Reusable Uretero-Renoscope




160K High Resolution



- ❖ Challenging Anatomy
- ❖ Passive Bending

Hugemed all in one system,
always be ready for portable image
processor



285°
Excellent
Deflection

For easier navigation
of lower pole



Flexible Sheath
3.6Fr working channel



Dual Buttons
For snap / Video capture



Self Locker
For reducing surgeon fatigue



True 1: 1
Torque ratio