



Video Gastroscope Video Colonoscopy

CMOS 1,000,000 Pixels HD resolution



Technical Specification

Item	Gastroscope	Colonoscope
Diameter of Distal-end	Ø9.8mm	Ø13.0mm
Diameter of Insert Tube	Ø9.6mm	Ø12.8mm
Biopsy channel	Ø 2.8mm	Ø 3.7mm
Working Length	1030mm	1330mm
Total Length	1330mm	1630mm
Depth of View	3-100mm	3-100mm
Resolution	1,000,000 CMOS	1,000,000 CMOS
Tip Deflection	U:210° D:900° L/R:100°	U:180° D:180° L/R:160°
Remark: We can provide OEM service , techical details can be customized.		

China Care Medical Equipment Co.,Ltd
 Adress:Room B101, Zhongbang Business Building, PuXing
 Road No. 87, Tianhe, Guangzhou,
 Chinaweb:www.chinacaremedical.com