



Handheld Digital Otoscope



Quality, Professional, Innovative

Handheld Digital Otoscope


Product Features

1. Super easy to operate. Focus and brightness are adjustable automatically, no need to manual adjustment.
2. A finger touch for taking photo or playing review, zoom in-out of 8 times, left/right movement.
3. With a built-in USB connection port, photos can be uploaded to computers at any time.
4. With a rechargeable Lithium-ion battery, ensure 4~5hrs' continuous working.
5. Special optical design allows a larger field of view (up to 40 degrees) to let doctors see what they want by eliminating the "searching" time by hooking the camera inner the ear.
6. 2.8" TFT color LCD with ultra-high resolution of 2560*1920.
7. Mobility of data storage; with a Micro SD Card still image (JPEG) and video (AVI) can be transferred to any computer.
8. Lens made of highest-grade miniature camera with 100,000 pixels.
9. Date and time displays on screen to record what you see in real time.
10. External LCD is optional which can be connected to otoscope to display larger images in real time.

Product Specifications

1. Intended use: Digital photograph for ear canal
2. Focus range: 20~60 mm
3. Dimension: 22*9*4 cm
4. Weight: 220 grams
5. Resolution: Full HD True 2560*1920
6. Light source: Natural white LED
7. White balance: Auto color calibration
8. LCD Display: 2.8" full color LCD
9. File format: JPEG (still picture), AVI (Video)
10. Interface: Mini USB, AV output
11. Memory: Micro SD card
12. Power: Li-ion rechargeable battery
13. Charge: 100V~240V AC 50/60Hz DC IN +5V 1A

Original Images



Product Appearance



China Care Medical Equipment Co.,Ltd
Address: Room B101, Zhongbang Business Building, PuXing
Road No. 87, Tianhe, Guangzhou,
China web: www.chinacaremedical.com