

VIDEO LARYNGOSCOPE



Features

1. Interchangeable Battery
2. One Key Control
3. Multi-Compatible Blades
4. Water Resistance
5. Direct-Power Cap

Specification

1. High resolution 3" LCD display monitor with 640×480 pixels for optimal view.
2. The display monitor is movable via 180° rotation axis.
3. Video imaging is in 4:3 format.
4. Rechargeable 18650 Li-ion Battery, 3.7 VDC, 3400mA.
5. With a fully charged 18650 battery, the VL's operating time should last 4+ hours.
6. Image device: Super Color VGA Sensor(640*480)
7. Lens View Angle (Field Of View): $55 \pm 10^\circ$



Configuration

Components	Qty/Unit	Model & Sizes	Remarks
Display monitor	1 piece	3"	
VL Handle	1 piece	M234	For Adults
18650 Li-ion Battery	2	3.7 VDC 3400 mAh	Rechargeable
Data cable	1	USB	
Charger Pak	1	For 18650 Use only	
Suitcase	1		Water Resistance
Direct-Power Cap	1		
Single Use VL Blade	5 pieces each size	Mac2/Mac3/Mac3D/MAC4	

Additional Note

1. The single use blade is made of imported medical grade PC from German.
2. In case of the electronic parts' failure, the displaying monitor does not work, by replacing it with the "Direct-Power Cap", the Video Laryngoscope can be used as a direct laryngoscope.
3. Water Resistance: the VL handle is IPX7, the Displaying Monitor is IPX2.
4. One Key controls all functions such as On/Off, Photo-taking and Video-Recording.
5. One VL handle is compatible with 4 different size blades, Mac2/Mac3/Mac4/Mac3d
6. With specially designed curve, Mac3D is suitable for difficult airway.

China Care Medical Equipment Co.,Ltd

Address: DaXin Building 2203C room, DeZheng North Road No.538, YueXiu District, Guangzhou, China. 510045

Cellphone/Whatsapp: 008613926480041 Skype: chinacaremed Wechat: Joseph20130708

Facebook/Youtube/Instragram/LinkedIn: china care medical

Website: www.chinacaremedical.com Email: chinacaremed@gmail.com