



4000A

— Hemodialysis Equipment —





The Power of Technology

Intelligent Operation

- 12.1 inch LCD+touch screen, multi-language operation interface.
- Adjust the fluid level of venous chamber to automatically prevent blood pollution and reduce dialysis risk.
- Convenient process for sequential dialysis (Dialysis $\leftrightarrow$ IUF), and high and low sodium(Na+) sequential dialysis.



Safe & Effective Treatment

- Mandatory self-test of the hydraulic, blood and monitoring systems to ensure the safety of treatment.
- Ultrasonic and optical detectors to ensure the safety of treatment.
- Dialysate feedback control system to ensure accurate dialysate concentration.
- Advanced balance feedback control system ensures UF precision.



Visual User Interface

- Preset automatic switch-off and disinfection mode, and record all disinfection activities.
- Dynamic tracking and monitoring of arterial pressure, venous pressure and transmembrane pressure to respond to abnormal situations quickly and trigger alarm timely.




Highly Compatible

- Operatable with various filter brands.
- Compatible with various brands of dialysis powder, bloodlines and other hemodialysis consumables.




Intelligent System Safe Treatment



Treatment modes

- HD
- Sequential HD
- IUF




RCM*

- Provides indicators such as urea clearance, and urea concentration decrease during the treatment.



BPM*

- Measures the blood pressure and heart rate of the patient, reducing workload of the medical staff.




Real-time Leak Monitoring

- Real - time leak monitoring system ensures the safety of treatment.



Back-up Power


- Blood circuit runs for more than 40 minutes.



* optional features


LED Alarm

- The LED alarm has three types of signals: red, yellow and green. Positioned on the top of the IV pole, visible from afar.



Multiple Profiles, Personalized Treatment

- Sodium profile
- UF rate profile
- Bicarbonate profile
- Heparin flow profile
- Dialysate flow profile
- Dialysate temperature profile
-



Bicart Holder

- Use Bicart online to avoid microbial contamination of B fluid.



Information System RJ45 Interface

- The information management system can be connected through the RJ45 interface to support refined management of the entire treatment process.



Expert Technical Solutions

- Our professional technical services team is dedicated to delivering prompt responses to hospital inquiries and needs. With the ability to offer solutions within 12-hours, they ensure a seamless user experience backed by a dependable technical support.
- Complimentary Ongoing Training: We provides regular training sessions at no cost for medical personnel, aimed at enhancing their understanding of machine operations and improving their responsiveness in emergency situations.





China Care Medical Equipment Co.,Ltd

Cellphone/Whatsapp: 008618998832732
Website: www.chinacaremedical.com

Wechat: Joseph20130708
Email: info@chinacaremedical.com