



**Skype:** chinacaremed

**Facebook:** China care medical

**Wechat:** joseph20130708

**Whatsapp:** +86 139 2648 0041

**Email:** joseph@chinacaremedical.com

**Website:** www.chinacaremedical.com

**Add:** Room B205, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China

---

## Biological Microscope

### Model: CCM-107B3

- Viewing Head: Sliding Trinocular Head with CCD adapter, Inclined at 45°,  
Interpupillary distance 55-75mm
- Eyepieces: Wide Field Eyepiece WF10X/18
- Objectives: Plan achromatic objectives: PL4X/NA0.1,  
PL10X/NA0.25, PL40X(s)/NA0.65,  
PL100X(s, oil)/NA1.25
- Nosepiece: Outward Quadruple Nosepiece
- Stage: Double layer Mechanical Stage ,Size:140X140mm /Movement range75X45mm
- Focus: Coaxial Coarse and fine Adjustment, Coars range 30mm ,  
fine range 0.002mm
- Condenser: Abbe N.A.=1.25 Condenser with iris diaphragm and blue filter
- Illumination: AC 220/110V,6V/20W Halogen lamp ,Brightness adjustable ,  
Plan-concave Mirror
- Package: 1set/1CTN,Foam+Carton
- Dimension: 32X27X43CM
- G.W.: 8KG/CTN,
- N.W.: 7KG/CTN

