



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

150L 250L Biochemical Incubators

Model: CCL-150A, CCL-250B



Introduction

This biochemical incubator is an indispensable thermostatic equipment for research and teaching in plant, biology, microbes, heredity, virus, iatrology and environment protection etc. with two-dimensional temperature controlling system. It is used widely in experiments of low and constant temperature, culture, and environment. The main characteristics are as follows:

1. The heat insulating material of incubator is made of foam frothed from polyurethane. It has strong ability of anti-interference with external hot or cool source.
2. The inner of incubator is made of engineering plastics by shaping techniques, so it has strong ability of corrosion resistance.
3. The working inner can be observed conveniently through the glass door.
4. To protect the refrigerating compressor, electrical outage and 4 minutes delay are designed in control line.
5. The temperature can be controlled automatically. Red LED display can show figure with clarity.

Technical Specification



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

Name	CCL-150A	CCL-250B
Specifications		
Volume	150 L	250 L
Range of temperature	5-50℃	5-50℃
Power of compressor	140W	150W
Heating power	200W	300W
Power supply	AC 220V 50Hz AC 110V 60HZ	AC 220V 50Hz AC 110V 60HZ

Note: Refrigerating when the ambient temperature at 16-32℃