




Auto-Cut extra edge



TAIWAN medical class pvc, more durable and high elastic, 12 month warranty



one-time forming mattress



Flaps fully wide



100% running testing 24 hours before shipment




China Care Medical
Equipment Co., Ltd.
Address: Room B101,
Zhongbang Business
Building, No. 87 Puxing
Road, Tianhe,
Guangzhou, China Web
: www.chinacaremedical.com



2.5'' Bubble Medical Air Mattress

Low noise pump with 18 month warranty
SGS approved
fire retardant Taiwan Nanya PVC



Specifications:

Pump	Pad
Item: CCFM012	Item: CCFM01
Power supply: 110V/220V,50-60Hz	Size: 200cm(L) x 90cm(W) x 7.0cm(H) (Standard Size) 190cm(L) X 82cm(W) X 7.0cm(H) (Econimic Size)
Air output: 5-7 liter/min*	130 bubble with end flaps
Pressure range: 50mmHg - 120 mmHg	Accessories: repair kit, 1.5M hose
Cycle time: 12 min*	Material: medical class Taiwan Nanya PVC K80
Power consumption: 7 W	PVC thickness: 0.3 mm
Fuse: 1A	Color: beige, blue, dark blue
Size: 23.5(L) x 12(W) x 9.5 (H)cm	Mattress weight: 1.8-2.1kg
Weight: 1.4 kg	Weight Capacity: 120-130kg