

Syringe Pump



TECHNICAL SPECIFICATIONS		
Features 4.3" Large LCD screen, with different colors display. - With LCD backlight, can use in various lighting conditions. - Time display, convenience for display and time memory. - Simultaneous display: time, battery, injection state, mode, speed, injection volume and time, syringe size, alarm sound, block, accuracy, body weight, drug dose and liquid amount. - Remote control function, can adjust speed, time, volume and drug amount directly, easier operation, save doctor and nurse's time. - Multi injection modes: Volume mode, Time mode, and Sequence mode. - Main interface: Three injection modes, Drug library, History record, Setting.	Purge Rate 50ml syringe: 1800ml/h (Bolus rate: 1200ml/h) 30ml syringe: 900ml/h (Bolus rate: 600ml/h) 20ml syringe: 600ml/h (Bolus rate: 400ml/h) 10ml syringe: 300ml/h (Bolus rate: 200ml/h) 5ml syringe: 150ml/h (Bolus rate: 100ml/h)	110~220V, 50/60 Hz, Built in rechargeable battery, 11.1V. Capacity >2000mAh. - After charging 8 hours, it can work at 5ml/h about 6 hours.
		Atmosphere Pressure 860hPa~1060hPa
		Working Surrounding - Temperature: +5℃ ~ +40℃ - Humidity: 20%~90%
		Storage Surrounding - Temperature: -30℃~+55℃ - Humidity: <95%
Technical Specification - Rate Range: 50ml syringe: 0.1~1800ml/h (0.1~999.9ml/h step: 0.1ml/h) 30ml syringe: 0.1~900ml/h (1000~1800ml/h step: 1ml/h) 20ml syringe: 0.1~600ml/h 10ml syringe: 0.1~300ml/h 5ml syringe: 0.1~ 150ml/h - Preset Volume Range: 0.1~2000.0ml (step: 0.1ml) - Maximum accumulated volume: 2000.0ml (step: 0.1ml) - Accuracy: Within ±2%	Alarm Injection finished, Preset dose injection finished, Occlusion, Syringe off, Syringe dislocated, Low battery, AC power off.	Safe Class Itype CF devices, internal battery, international load continuous working.
	Occlusion Three levels for acclusion: High, Middle, Low.	Size Size: 282mm(L)x210mm(W)x148mm(H)
	KVO Rate 0.1~5ml/h (step: 0.1ml/h)	Weight Weight: 2.1kg
	Fuse F1AL/250V("F" indicate fast, low breaking capacity), 2pcs (installed inside pump)	
	Power Supply	

NOTE: The specifications are subject to changes without prior notice.



China Care Medical Equipment Co., Ltd.
Address: Room B101, Zhongbang Business Building, No. 87 Puxing Road, Tianhe, Guangzhou,
China Web: www.chinacaremedical.com

- 4.3" color LCD screen, backlight display, can be used in various lighting conditions.
- Advanced technology, based on Linux system, more safe and steady.
- Visible and audible alarms cover all abnormal conditions.
- Drug library, History record, Remote control.



Syringe Pump

Screen

- 4.3" Color TFT LCD screen.
- With high resolution with LED backlight.



Syringe

- Automatic identification function of syringe specifications:
It can be used for syringes of 5ml, 10ml, 20ml, 30ml and 50/60ml.




Remote Control


- More convenient and save time.
- Remote control is instead of button operation.
- The user can control the infusion pump within 5 meters.




Syringe Pump



Volume Mode



Time Mode



Sequence Mode

Multi injection modes

- Volume mode(Rate mode);
- Time mode;
- Sequence Mode.



The back
The power port has a dust-proof cover

