



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China

Biological Binocular Microscope



Model: CCM-207

Viewing Head:	Compensation Binocular Head, Inclined at 45°, Interpupillary distance 55-75mm
Eyepieces	Wide Field Eyepiece WF10X 16W
Objectives	Achromatic objectives 4X 10X 40X(S) 100X(S,Oil)
Nosepiece	Outward Quadruple Nosepiece
Stage	Double layer Stage, Size: 140X140mm
Focus	Coaxial Coarse and fine Adjustment, Coars adjustment 30mm ,fine range 0.002mm
Condenser	Abbe N.A.=1.25 Condenser with iris diaphragm and filter
Illumination	AC 220/110V, 6V/20W Halogen lamp ,Brightness adjustable ,plan-concave mirror
Package:	1set/1CTN, foam+carton+wooden case
Dimension:	38X33X53CM
G.W.:	12KG/CTN,
N.W.:	8KG/CTN