



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

Color LCD Infusion Pump

Model: CCS-1200





Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China



Features:

1. Compatible with IV sets of any brand, universal type.
2. Automatically calibration for IV sets of any brand.
3. Three working modes: Drop counting control, Rate control and Time control mode.
4. 6 Hours continuous work when the battery is fully charged.
5. Have a unique color LCD and an internal battery with a high capability.
6. Switch between ml/h and drops/min freely.
7. Three levels of alarm volume: high, medium and low.
8. Three levels of occlusion pressure: high, medium and low.
9. KVO function: When finish the total volume of an infusion, the pump will alarm and enter KVO state.

Technical Specifications

PARAMETER	SPECIFICATION
Drop Detector	Outside caliper rotary detector (patent) Infrared sensor



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China

Control Way	Drop counting and adjusting the delivery rate automatically
Fit IV tubes	20d/ml,60d/ml and other, as well as brown tubes
Pump Type	Peristaltic finger pump
Accuracy	±2%, when delivery rate ≤25drop/min (drop counting & time control mode) ±3%, when delivery rate >25drop/min (drop counting & time control mode) ±5% (Volume control mode after the correct calibration)
Delivery rate range	1-1999min (1min step)
Volume limit range:	0-9999ml (1 ml step)
Total Volume Range	0-9999.9ml (1 ml step)
Purge Rate	600.0ml/h
KVO Rate	KVO is 1ml/h when the delivery rate≤10ml/h, Otherwise KVO is 3ml/h.
Occlusion pressure	Three level setting LOW, MEDIUM or HIGH
Alarms	Near empty, Syringe empty, Volume complete (when limit volume is set), Occlusion, Syringe loose, error, Low battery, device failures. Three levels of alarm volume LOW, MEDIUM or HIGH
Bolus function	Bolus rate:10-1200ml/h,Bolus volume: maximum of 100ml
Power source	AC power supply: AC 100 to 240V(+/-10%), 50/60Hz; Built-in Lithium Battery: Rechargeable,11.1V,2000mAh. 6 hours of continuous operation
Accessories:	AC power line, Operational manual
Protection Class	Class I, Type CF,IPX4
Dimension /Net weight	130mm*200mm*230mm , 2.5kg