



Healthy Every Day

China Care Medical Equipment Co.,Ltd

Room B101, ZhongBang Business Building,
PuXing Road No. 87, Tianhe, Guangzhou, China

+86 139 2648 0041


info@chinacaremedical.com



HealFusion

Syringe Pump s5/s5-plus
Infusion Pump i5/i5-plus





It applies sophisticated stepping motor and driving system and the stepping rate is 0.01ml/h, which guarantees the accuracy; Faulty tree analysis optimizes every step and assures the stability and safety of our pumps. Advanced touch screen technology makes all operations convenient and intelligent that can help to ease the workload of nurses. Moreover, with its small volume and lightweight, abundant infusion modes meet various clinic demands and improve its benefit without complicated conversion and operation.

>>> Convenient

- ◆ Modular design, abundant infusion modes meet various clinic demands without complicated conversion and operation.
- ◆ It is convenient for transportation of patients when it is installed on the elaborate handle.
- ◆ The shell is made from PBT plastic and it has well characteristics of resistance to corrosion, which is more convenient to clean and disinfect.

>>> Easy to operate

- ◆ It applies large touch screen and it is easy to operate, which is more convenient and time-saving than the traditional keyboard operation.
- ◆ Infusion mode: rate mode, time rate mode, preset time mode, Dose mode (not includes i5/s5); all of these meet various clinic demands without complicated conversion and operation.

>>> Accurate

- ◆ It applies sophisticated stepping motor and driving system and the stepping rate is 0.01ml/h, that can totally meet clinic demands of toxic and anesthetic infusion or the high accuracy neonatal infusion.
- ◆ Syringe pump has fast start functions, which can eliminate error caused by starting delay.
- ◆ Infusion pump has velocity compensation function, which can avoid the error caused by pipeline collapse.

>>> Safe

- ◆ Dual CPU design assures the safety of system.
- ◆ Online titration function can directly modify rate without stopping.
- ◆ Anti-bolus function avoids the adverse reaction caused by blood concentration.
- ◆ Dynamic Pressure Surveillance (DPS): 5 levels of alarm range is adjustable, it not only helps medical staff to set reasonable pressure alarm range, but also is helpful to find out pipeline off, leakage of injection and clogging in time.

