



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joechinacare@gmail.com

Website: www.chinacaremedical.com

Add: DaXin Building 2203C room,DeZheng North Road No.538,YueXiu District,Guangzhou,China



Viewing Head: Sliding Binocular Head, Inclined at 45° ,
Interpupillary distance 55-75mm
Eyepieces Wide Field Eyepiece WF10X 16W
Objectives Achromatic objectives 4X 10X 40X(S) 100X(S,Oil)
Nosepiece Outward Quadruple Nosepiece
Stage Double layer Stage, Size: 140X140mm
Focus Coaxial Coarse and fine Adjustment,
Coars adjustment 30mm , fine range 0.002mm
Condenser Abbe N.A.=1.25 Condenser with iris diaphragm and filter
Illumination AC 220/110V, 6V/20W Halogen lamp , Brightness adjustable , plan-concave mirror

Package: 1set/1CTN, foam+carton+wooden case
Dimension: 38X33X53CM
G.W.: 12KG/CTN,
N.W.: 8KG/CTN