



**Skype:** chinacaremed

**Facebook:** China care medical

**Wechat:** joseph20130708

**Whatsapp:** +86 139 2648 0041

**Email:** joseph@chinacaremedical.com

**Website:** www.chinacaremedical.com

**Add:** Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

---

## Microplate Mixers

Model: **CCL-MXM**



### Features

- 1, Wide range of accessories
- 2, A wide range of speed from 0 to 1500rpm
- 3, Used for various mixing applications with optional adapters
- 4, Maintenance-free brushless DC motor

### Specifications

Shaking movement: Orbital

Orbital diameter: 4.5mm

Max. shaking weight (with accessory): 0.5kg

Motor type: Brushless DC motor

Motor rating input: 18W

Motor rating output: 10W

Speed range: Single microplate: 0-1500rpm, Double microplate: 0-1000rpm

Speed display: Scale

Voltage: 100-240V, 50/60Hz

Power: 20W

Dimension: 260×150×80mm

Weight: 3kg

Permissible ambient temperature and humidity: 5-40 degree, 80%RH

Protection class: IP21