

CCL-MXS

Mélangeur Vortex à vitesse variable

Caractéristiques

- Vitesse variable de 0 à 3 000 tr/min pour une agitation douce ou un mélange rapide des échantillons
- Mode de fonctionnement tactile ou continu
- Comprend une tête supérieure standard pour mélanger un seul tube à essai ou un autre petit récipient jusqu'à 30 mm de diamètre
- Le contrôle électronique de la vitesse maintient une vitesse constante même en cas de changement de charge
- La base en acier avec des pieds d'aspiration spécialement conçus assure la stabilité et amortit les vibrations
- Construction robuste en fonte d'aluminium

VT.1.1 Standard top (par défaut)

Réf. 18900034

for <Ø30mm tubes and small vessels



Caractéristiques	CCL-MXS (Vitesse réglable)
Tension, fréquence	100-120V/200-240V,50/60Hz
Puissance	60W
Mouvement de mélange	Orbital
Diamètre orbital	4mm
Type de moteur	Shaded pole motor
Puissance nominale du moteur en entrée	58W
Puissance nominale du moteur en sortie	10W
Plage de vitesse	0-3000rpm
Affichage de la vitesse	Scale
Type de fonctionnement	Touch operation/Continuous
Dimensions (L × H × P)	127 × 130 × 160mm
Poids	2.0kg
Température ambiante admissible	5-40°C
Humidité relative admissible	80%RH
Classe de protection	IP20

China Care Medical Equipment Co.,Ltd

Room B101, ZhongBang Business Building, PuXing Road No. 87, Tianhe, Guangzhou, China

www.chinacaremedical.com



VT1.3.1 tube adapter	Cat. No.18900020	VT1.3.10 tube adapter	Cat. No.18900855
for 48 test tubes, Ø6mm (1mL)		for 4 holes test tubes, Ø28mm (50mL)	
VT1.3.2 tube adapter	Cat. No.18900021	VT1.3.11 Tube adapter	Cat. No.18900862
for 15 test tubes, Ø10mm (1.5-2mL)		12XØ9.8mm+12XØ6.6mm+12XØ5mm (1.5/2.0/0.5/0.2mL)	
VT1.3.3 tube adapter	Cat. No.18900022	VT1.2 Tube holding rod	Cat. No.18900044
for 16 test tubes, Ø12mm (5mL)		used with tube adapters	
VT1.3.4 tube adapter	Cat. No.18900023	VT1.3 Universal top plate	Cat. No.18900505
for 8 test tubes, Ø16mm (10mL)		Ø 100mm Includes 4 vacuum chuck	
VT1.3.5 tube adapter	Cat. No.18900024	PS1.0 Microplate clamp	Cat. No.18204342
for 8 test tubes, Ø20mm (10-15mL)		Includes rectangle PCR tube adapter	
VT1.3.6 Platform pad	Cat. No.18900043	VT1.3.12 Large Platform Multi-Specimen Test Tube Adapter	Cat. No.18900875
Ø99mm, multiple small tubes and vessels		44 holesXØ5mm+55 holesXØ6.6mm+30 holesXØ9.8mm (1.5/2.0/0.5/0.2mL)	
VT1.3.7 vacuum chuck	Cat. No.18900158	VT1.4 Ampule/Tube Attachment	Cat. No.18204566
made of rubber		4pc ampule/tube (0.2mL-15mL)	
VT1.3.8 tube adapter	Cat. No.18900395	VT1.5 Erlenmeyer flask attachment	Cat. No.18204561
for 24 holes test tubes, Ø 10mm (1.5-2mL)		100mL-250mL, Operating speed: 0~1500rpm	
VT1.3.9 PCR tube adapter	Cat. No.18204185		
for 56 test PCR tubes, strips of 8, Ø5mm (0.1mL-0.5mL)			





Application d'accessoires

Accessories	Modèle à vitesse réglable (0-2500rpm)	
	Mode tactile (zone haute vitesse)	Mode continu (zone basse vitesse)
VTI.1	Y	Y
VTI.2 + VT1.3.1	Y	
VTI.2 + VT1.3.2	Y	
VTI.2 + VT1.3.3	Y	
VTI.2 + VT1.3.4	Y	
VTI.2 + VT1.3.5	Y	
VTI.3 + VT1.3.6	Y	
VTI.3 + VT1.3.1		Y
VTI.3 + VT1.3.2		Y
VTI.3 + VT1.3.3		Y
VTI.3 + VT1.3.4		Y
VTI.3 + VT1.3.5		Y