



China Care Medical Equipment Co., Ltd.

Address: Room B101, Zhongbang Business Building, No. 87 Puxing Road, Tianhe, Guangzhou, China Web: www.chinacaremedical.com



Vortex Mixer is designed for mixing liquids in schools, laboratories and factories. It is not suitable for use in residential area or other areas that may cause danger to the user or instrument. An ideal experiment tool for laboratory in schools, environment protection, research, sanitation, immunisation, petrochemical, metallurgy, chemical and foodstuff organizations.

CCL-XH-C

Technical data:

power supply: AC220V 50–60HZ

shaking frequency: 0–2800rpm

stepless regulation

Power: 60W

rubber head, continuous

Shake mode: 