



Veterinary portable dental X-ray imaging solution

- USB Sensor Port
- Long Battery Life
- Wireless Charging
- All-in-One Design
- Instant dental imaging
- Built-in imaging system
- Low-dose portable X-ray
- Nano lead shielding glass






Veterinary portable dental X-ray imaging solution



Type	Gun-type
Special function	Instant dental imaging
Tube Voltage	70.0KV
Exposure time range	0.01-2.0S, adjustable
Battery	DC22.0V/700.0mAh
Invert Frequency	200.0kHz
System	Human/Veterinary,adjustable

Screen	2.8' touch screen
Special configuration	Wireless charging case
Tube current	3.5mA
Focal Spot to Skin Distance	≥200.0mm
Voltage	110.0-220.0V,50.0-60.0Hz
Compatible Sensors	iRay, Brigates , Handy,etc.
Weight	1.8kgs

High-definition dental images



Standard and Optional Features



China Care Medical Equipment Co., Ltd.
Address: Room B101, Zhongbang Business Building, No. 87 Puxing
Road, Tianhe, Guangzhou, China Web: www.chinacaremedical.com