



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China



Viewing Head:	Compensation Binocular Head,Inclined at 30°, Interpupillary distance 55-75mm
Eyepieces	Wide Field Eyepiece WF10X
Objectives	Achromatic objectives 4X 10X 40X(S) 100X(S,Oil)
Nosepiece	Outward Quadruple Nosepiece
Stage	Single Stage with attachable ruler ,Size:115X125mm
Focus	separate coarse and fine adjustment, coarse adjustment focusing15mm. Fine adjustment 1.6mm
Condenser	Abbe N.A.=1.25 Condenser with iris diaphragm and filter the center of condenser is adjustable
Illumination	AC 220/110V,6V/20W Halogen lamp ,Brightness adjustable ,plan-concave mirror
Package:	1set/1CTN,foam+carton
Dimension:	41.5X31.5X21.5
G.W.:	6KG/CTN,
N.W.:	5KG/CTN