


Veterinary Dental X-Ray Machine



- High level of safety protection
- Operation with Remote controller
- Ultra-high frequency X-ray generator
- Top tube voltage and tube current parameter values
- Tube voltage and tube current can be manually adjusted
- The world's leading ultra-high frequency inverter technology






Veterinary Dental X-Ray Machine




Type	Gun-type
Tube Voltage	60.0,65.0,70.0KV (adjustable in three levels)
Tube current	2.5,3.5mA(adjustable in two levels)
Exposure time	0.01- 1.0S
Tube focus	0.4mm
Distance from skin to cone	200.0mm
System	Human, Veterinary
Language	English, Chinese, Russian
Launch tube	Short tube, Long tube

Screen	2.8' LCD touch screen
KV Ripple Frequency	<1.0%
Inverting Frequency	>100.0KHZ
Radiation dose	2.447mGy
Surface leakage radiation dose	0.0007mSv
Lithium battery capacity	DC22.2V/950.0mAh
Power Supply	AC100.0-240.0V,50.0/60.0HZ
Net Weight	2.0KGS
Gross Weight	4.0KGS
Carton Size	41.00*34.50*18.7cm

Excellent imaging



Related products



China Care Medical Equipment Co., Ltd.
Address: Room B101, Zhongbang Business Building, No. 87 Puxing Road,
Tianhe, Guangzhou,
China Web: www.chinacaremedical.com