



Dental X-Ray Machine

- Ultra high frequency X-ray generator
- 1.0kg class with reducing weight by 50.0%
- Ultra-low radiation dose
- High level of safety protection
- The world's leading ultra-high frequency inverter technology







Dental X-Ray Machine




Type	Camera-type
Screen	LCD Touch screen
Tube Voltage	60.0KV, 65.0KV, 70.0KV adjustable
Tube current	2.5mA, 3.5mA adjustable
Tube focus	0.8mm
Focus-skin distance	200.0mm
Exposure time	0.010-1.0S

KV Ripple Frequency	< 1.0%
Inverting Frequency	> 100.0KHZ
Input	100.0-240.0V AC, 50.0-60.0HZ
Lithium battery capacity	DC22.2V/950.0mAh
Net Weight	1.8KGS
Gross Weight	3.0KGS
Carton Size	34.0cm*34.0cm*15.0cm

Excellent imaging



Related products



China Care Medical Equipment Co., Ltd.
Address: Room B101, Zhongbang Business Building, No. 87
Puxing Road, Tianhe, Guangzhou, China Web: www.chinacaremedical.com