



## Dental X-Ray Machine

- Ultra high frequency X-ray generator
- 1.0kg class with reducing weight by 50.0%
- Ultra-low radiation dose
- High level of safety protection
- The world's leading ultra-high frequency inverter technology







## Dental X-Ray Machine




|                     |                                   |
|---------------------|-----------------------------------|
| Type                | Camera-type                       |
| Screen              | LCD Touch screen                  |
| Tube Voltage        | 60.0KV, 65.0KV, 70.0KV adjustable |
| Tube current        | 2.5mA, 3.5mA adjustable           |
| Tube focus          | 0.8mm                             |
| Focus-skin distance | 200.0mm                           |
| Exposure time       | 0.010-1.0S                        |

|                          |                              |
|--------------------------|------------------------------|
| KV Ripple Frequency      | < 1.0%                       |
| Inverting Frequency      | > 100.0KHZ                   |
| Input                    | 100.0-240.0V AC, 50.0-60.0HZ |
| Lithium battery capacity | DC22.2V/950.0mAh             |
| Net Weight               | 1.8KGS                       |
| Gross Weight             | 3.0KGS                       |
| Carton Size              | 34.0cm*34.0cm*15.0cm         |

## Excellent imaging



## Related products



China Care Medical Equipment Co., Ltd.

Address: Room B101, Zhongbang Business Building, No. 87 Puxing Road, Tianhe, Guangzhou,  
China Web: [www.chinacaremedical.com](http://www.chinacaremedical.com)