



Vet

Veterinary Dental X-Ray Machine

- Ultra high frequency X-ray generator
- 1.0kg class with reducing weight by 50.0%
- Ultra-low radiation dose
- High level of safety protection
- The world's leading ultra-high frequency inverter technology





Vet

Veterinary Dental X-Ray Machine




Type	Camera-type	KV Ripple Frequency	< 1.0%
Screen	LCD Touch screen	Inverting Frequency	> 100.0KHZ
Tube Voltage	60.0KV, 65.0KV, 70.0KV adjustable	Input	100.0-240.0V AC, 50.0-60.0HZ
Tube current	2.5mA, 3.5mA adjustable	Lithium battery capacity	DC22.2V/950.0mAh
Tube focus	0.8mm	Net Weight	1.8KGS
Focus-skin distance	200.0mm	Gross Weight	3.0KGS
Exposure time	0.010-1.0S	Carton Size	34.0cm*34.0cm*15.0cm

Excellent imaging



Related products



China Care Medical Equipment Co., Ltd.
Address: Room B101, Zhongbang Business Building,
No. 87 Puxing Road, Tianhe, Guangzhou, China Web:
www.chinacaremedical.com