

Ophthalmic Wide Field Imaging System

CC25041709

Portable compact device capturing high-resolution images and the benefit of video recording. Simple to operate and a USB 3.

0 compatible device can easily work with any laptop. Suitable for various screening scenarios.

Technical data

CC25041709 is a WideField camera, very light with easy handling. Is integrating Light Source, Refractive compensation and High resolution Imaging system

Probe Diameter	8.1mm (smallest probe in the market)
Probe Weight	527.3gr (the lightest in the market)
Probe connectivity	USB3.0 Plug & Play
Wide Field of View	135° to 144°
Light Source	Color Light Source; Fluorescence Angiography
Center field of view	≥ 30 lp/mm
Middle field of view	$(\pm 22.5^\circ) \geq 20$ lp/mm
Edge field of view	$(\pm 45^\circ) \geq 15$ lp/mm
Illumination Source	white LED