



Pharmacy Refrigerator



- ▶ The pharmacy refrigerator is the specific medicine cold-storage equipment. It can be also used for storing reagent, biological product, etc. It is especially available for pharmacy, pharmaceutical factory, epidemic station, blood bank station, health center and lab.

Pharmacy Refrigerator (368L)

Specifications

Model	
Temperature Range (°C)	2~8°C
Capacity (L/Cu. Ft)	368L
Chamber	6
Cooling Method	Forced Air Cooling System
Defrost	Auto
Refrigerant	R600a CFC Free
Power consumption	2.5KW.H/24H
Ambient Temperature (°C)	10-32
Temperature Controller	Microprocessor
Display	LED Digital Display
Inside Material	Color Sprayed Steel
Outside Material	Color Sprayed Steel
Outside Door/Type	Single Glass Door
Insulation	PURF
Alarm	Power failure alarm function High/Low Temperature, Power Failure, Sensor Failure, Door Ajar, System Failure Alarm
Alarm Type	Sound/Light
Castor	Yes
Lamp	LED
Noise (db)	42
Voltage/Frequency (V/Hz)	220V/50Hz
Input Power (W)	250
Climate type	N
Electric shock resistance	Class 1 Type B
Exterior (DxWxH) (mm)	600*600*1940mm
Package Dimension (DxWxH) (mm)	730*730*2060mm
Net Weight (kg)	75
Gross Weight (kg)	85
Optional function	Portable temperature recorder with USB interface Temperature record printer Temperature record printer Remote alarm (wireless) UPS emergency power Solar energy power source DC24V-AC220V inverter



China Care Medical Equipment Co., Ltd.
Address: Room B101, Zhongbang Business Building, No.
87 Puxing Road, Tianhe, Guangzhou, China Web: www.chinacaremedical.com