


# CC25042905

## Nd:YAG Laser for Ophthalmology



## Specifications

- Laser Type: Nd:YAG, Q-Switched
- Wavelength: 1064 nm
- Burst: 1,2,3 pulses per burst
- Pulse Duration: 4.5 ns
- Repetition Rate: 2.5Hz
- Energy (Power)
  - A) 1 pulse maximum output: 11mJ
  - B) 2 pulse maximum output: 19 mJ
  - C) 3 pulse maximum output: 28 mJ
  - D) Output Energy Regulation: 100% ~ 6% adjustable by 7 steps
- Cone Angle: 18°
- Spot Size: 30  $\mu\text{m}$
- Aiming Beam Wavelength: 633 nm
- Cooling System: Air Cooled
- Dimension: 765mm x 610mm x 420mm
- Weight: 19Kg



### China Care Medical Equipment Co.,Ltd

- 📍 Room B101, ZhongBang Business Building,  
PuXing Road No. 87, Tianhe, Guangzhou, China
- ☎ +86 189 9883 2732
- ✉ info@chinacaremedical.com
- 🌐 www.chinacaremedical.com