



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: info@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room C208, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China



Technical Indicators

operating wavelength: 1064nm

pulse energy: 0-40mJ

width of pulse: 6 ns



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: info@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room C208, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

focal aperture angle: 16 angle

repetition frequency: 2.5 Hz

number of pulses: 1-4 pulse can be randomly taken

aiming beam: double red laser spot

energy selection: 100%、70%、50%、35%、25%、12%、6%

cooling-down method: cooling-down by air

power requirement: 220V/50Hz、500W