



Ophthalmoscope(DC)
Model:CCE- YZ11

The product is elegant and portable.

Detailed profile

Characteristics:

- Elegant appearances.
- It is convenient for mobile inspection.Both common dry cell and rechargeable battery power are applicable.

Specifications:

Illumination Form	Large Spot, Middle Spot, Small Spot, Red Free
Diopter Compensation	0 ±1 ±2 ±3 ±4 ±5 ±6 ±8 ±10 ±12 ±16 ±20 -25 -35D
Illumination Source	2.6V/0.78W,Incandescent Bulb
Power	2#Alkalescent battery ,rechargeable battery