

Ophthalmoscope
Model:CCE- YZ6E



Observed images of the fundus are sharp and upward

Characteristics:

- Observed images of the fundus are sharp and upward.
- Halogen bulb with bright and even illumination.
- Automatic turning off effectively protects the bulb.
- Two-step adjustable illumination: S(Bright),W(Dark).

Specifications:

Illumination Form	Large Spot, Small Spot, Slit, Net, Red Free
Diaoptor Compensation	-35D~+20D, 24 steps in total
Illumination Source	6V/5W, Halogen Bulb
Power Supply	AC220V/ 50Hz, AC110V / 60Hz

China Care Medical Equipment Co., Ltd.
Address: Room B101, Zhongbang Business Building, No. 87 Puxing Road, Tianhe
, Guangzhou, China Web: www.chinacaremedical.com