



Tonometer

Model:CCE- YZ7A

The product is designed and manufactured conforming to the international standard ISO/TR 8612---1997,which is reliable enough in quality,smart and portable in structure and convenient for use.

Detailed profile

The product is designed and manufactured conforming to the international standard ISO/TR 8612---1997,which is reliable enough in quality,smart and portable in structure and convenient for use.

Principle of design: use the pointer to press cornea and then measure the pressed depth ,thus indirectly get the intraocular pressure.

Characteristics:

- Schiotz type tonometer.
- All the parts contacting with cornea are made of stainless steel which is anti-acid and does no harm to issues.
- Safe, reliable, practical, accurate and portable.

China Care Medical Equipment Co., Ltd.

Address: Room B101, Zhongbang Business Building, No. 87 Puxing Road,
Tianhe, Guangzhou, China Web: www.chinacaremedical.com