

Ophthalmoscope retinoscope(LED Bulb)

Model:CCE- YZ6H



The product adopts the imported advanced LED bulb which is durable enough(with designed service life of 100,000 hours) and has even light spot.

Detailed profile

Characteristics:

- LED cold light source. It is durable enough with the designed service life of 100000 hours.
- Low heat benefits from the cold light source.
- Imported bulb maintains 3200K color temperature for long time.The illumination is soft and more similar with sunlight.
- Ultraviolet radiation is eliminated to fully protect patient's eyes.
- Even light spot and sharp observed fundus image.
- Soft rubber around the eyeholes protect eyes and stabilize the ophthalmoscope.

Specifications:

Illumination Form	Large Spot, Small Spot, Slit, Net, Red Free
Diopter Compensation	0 $\pm$ 1 $\pm$ 2 $\pm$ 3 $\pm$ 4 $\pm$ 5 $\pm$ 6 $\pm$ 8 $\pm$ 10 $\pm$ 12 $\pm$ 16 $\pm$ 20 -25 -35D
Illumination Source	3.5V/1W LED Bulb
Input Voltage	AC220V $\pm$ 22V/50Hz $\pm$ 1Hz,AC110V $\pm$ 11V/60Hz $\pm$ 1Hz

China Care Medical Equipment Co., Ltd.
Address: Room B101, Zhongbang Business Building, No. 87 Puxing
Road, Tianhe, Guangzhou,China Web: www.chinacaremedical.com