
Tonometer



Model:CCE- YZ30X

The product is designed on the principle basis of Goldman tonometer.It can be connected with slit lamp(Carl Zeiss type).

Detailed profile

Characteristics:

- It's not only used to examine the eyes,but also measure the intraocular pressure while working under the slit lamp.
- Accurate measurement ensures the total tolerance less than $\pm 0.066\text{KPa}$ (0.5mmHg).
- Directly get the intraocular pressure.
- Measured intraocular pressure will not be affected by hardness of eye.
- Adjustable measuring pressure ensures the working stability and accurate measurement of tonometer in long time.

China Care Medical Equipment Co., Ltd.

Address: Room B101, Zhongbang Business Building, No. 87 Puxing Road, Tianhe, Guangzhou, China Web: www.chinacaremedical.com