

Ophthalmoscope
Model:CCE- YZ11D



The portable product has excellent optical performance by improved optical design.

Detailed profile

Characteristics:

- Improved light path and excellent optical performance.
- Imported halogen bulb provides high and adjustable brightness and good color rendering.
- Smart, portable and convenient for mobile inspection.
- Equipped with special charger, quick and convenient recharging, against over recharging.
- Anti-recharging and anti-discharging function protects the battery.

Specifications:

Illumination Form	Large Spot, Small Spot, Slit, Net, Red Free
Diopter Compensation	0 ±1 ±2 ±3 ±4 ±5 ±6 ±8 ±10 ±12 ±16 ±20 -25 -35D
Illumination Source	3.5V/2.8W,Halogen Bulb
Input Voltage	AC220V±22V/50Hz±1Hz,AC110V±11V/60Hz±1Hz
Power(Main body of Ophthalmoscope)	Li+ Battery Pack 3.7V,Rechargeable