



China Care Medical Equipment Co., Limited

Tel: +86 20 62871996

Mob: +86 13711403953

Contact person: Shally Xie

Email: shallychinacare@gmail.com

Web: <http://www.chinacaremedical.com>

Add: DaXin Building 2203C room, DeZheng North Road No.538, YueXiu District, Guangzhou, China

Slit Lamp

Model: CCE-YZ5E



YZ5E with excellent performance and high quality is found very popular in the market.

Detailed profile

Characteristics:

- Advanced converging stereoscopic microscope.
- Clearer optical image and broader field of view.
- Optical lens are all moistureproof, mildewproof and anti-reflection.
- Equipped with Trihedron lens ,Gonioscope,Funduscope.
- Three-step magnifications for eyepiece are optional:10x,16x,25.5x.

For free: Trihedron lens ,Gonioscope,Funduscope

YZ5E is a basic type.

YZ5EIII is equipped with double-teaching tube.

Specifications:

Type	Converging stereoscope
Magnifications	Two-step magnifications
Eyepiece	10x,16x
Magnification Ratio(Field of View)	10x(D18),16x(D14.5) 16x(D11.25),25.6x(D9)
Range of Pupil Distance	10x:55mm~82mm 16x:51mm~78mm
Adjustment of Diopter	10x: $\pm 8D$ 16x: $\pm 10D$
Slit Width(mm)	0~9 adjustable(slit is round when slit width is 9mm)
Slit Height(mm)	1~8 adjustable

**China Care Medical Equipment Co., Limited****Tel:** +86 20 62871996**Mob:** +86 13711403953**Contact person:** Shally Xie**Email:** shallychinacare@gmail.com**Web:** <http://www.chinacaremedical.com>**Add:** DaXin Building 2203C room, DeZheng North Road No.538, YueXiu District, Guangzhou, China

Slit Angle	0-180° adjustable
Slit Inclination	0-23°
Aperture Diameter(mm)	D9,D8,D5,D3,D2,D0.2
Filter	Heat absorption, grey, redfree, cobalt blue
Illumination Bulb	12V/50W halogen bulb
Optional Accessories	Applanation Tonometer(T type)