



China Care Medical Equipment Co., Limited

Tel: +86 20 62871996

Mob: +86 13711403953

Contact person: Shally Xie

Email: shallychinacare@gmail.com

Web: <http://www.chinacaremedical.com>

Add: DaXin Building 2203C room, DeZheng North Road No.538, YueXiu District, Guangzhou, China

Slit Lamp Microscope

Model: CCE-YZ5G



The product with complete function and advanced design is well received in domestic and abroad market.

Characteristics:

- Professional digital camera is more sensitive to the changes of patient's eyes, even the smallest changes.
- 14mm aperture provides broader and clearer field of view.
- Imported Schott optics from Germany reduces optical distortion, increases clarity and depth of focus.
- All the Optical lenses are moistureproof, mildewproof and anti-reflected treated.
 - Osram bulb provides the low heating, high brightness and even lighting.
- Five-step drum magnification for eyepiece. Sclera vein could be observed under the 40x magnification
- Excellent optical performances and long span usage.
- Anti-ultraviolet radiation filter ensures patient's safety during examination.
- All the joints can be connected by accessories with international standards.

**China Care Medical Equipment Co., Limited****Tel:** +86 20 62871996**Mob:** +86 13711403953**Contact person:** Shally Xie**Email:** shallychinacare@gmail.com**Web:** <http://www.chinacaremedical.com>**Add:** DaXin Building 2203C room, DeZheng North Road No.538, YueXiu District, Guangzhou, China**Specifications:**

Microscope Type	Galilean stereoscopic microscope
Objective Magnification	Five-step Drum Magnification
Eye Piece	12.5x
Magnifications (Field of View)	6x (ϕ 33mm) 、 10x (ϕ 22.5mm) 、 16x (ϕ 14mm) 、 25x (ϕ 8.8mm) 、 40x (ϕ 5.5mm)
Pupil Distance	55mm~75mm
Diopter Adjustment	-5D~+3D
Slit Width	0mm~14mm adjustable(slit is round when the slit width is 14mm)
Slit Height	1mm~14mm adjustable
Slit Angle	0°~180° adjustable
Slit Inclination	5°、 10°、 15°、 20°
Diameter of Light Spot	ϕ 14mm、 ϕ 10mm、 ϕ 5mm、 ϕ 2mm、 ϕ 1mm、 ϕ 0.2mm
Filter	Heat Absorption, Grey, Red Free, Cobalt Blue
Input Voltage	AC220V $\pm$ 22 V; AC110V $\pm$ 11 V
Electrical Safety Standard	Executive Standard: IEC601—1, Class I, Type B
Packing Volume	720mmx495mmx480mm
Total Weight	24Kg
Optional Accessories	Applanation Tonometer 、 Measuring Eye Piece , Trihedron Lens、 Gonioscope、 Funduscope、 Aspheric Lens、 Teaching Tube