

Keratometer

Model:CCE-YZ38



The product is applied in measuring diopter on the surface of cornea, and the axis and diopter of cornea astigmatism. Also, it can measure the diopter of implanting lens assisted by A ultrasonic.

Detailed profile

Specifications:

Measuring Range	Curvature Radius	5.5—11mm
	Diopter	30—60D
Minimum Division	Radius	0.02mm
	Diopter	0.25D
Minimum Surface Area Required by Measurement	r=5.5mm	D1.65mm
	r=7.5mm	D2.36mm
	r=11mm	D3.36mm

China Care Medical Equipment Co., Ltd.

Address: Room B101, Zhongbang Business Building, No. 87

Puxing Road, Tianhe, Guangzhou, China Web: www.chinacaremedical.com