



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

PORTABLE HEALTHCARE OXYGEN CONCENTRATOR



Model: CCS-JAY1

Feature

1. Continuous flow, light weight.
2. It can be run by car adapter DC12V, Lithium battery (one hour) 220V-240V or 110V

Application

Especially for shopping, traveling, etc. outdoor application

Patients' short time outdoor use.

Headache dizziness, sleep apnea caused by cervical spondylosis

Mental worker and student recover energy

Pregnant and the aged to keep healthy

Oxygen cocktail making

Specifications

Model: CCS-JAY1

Flow Rate: 1-5LPM

Purity: 93%~38%

Outlet Pressure: 0.03-0.07MPA

Sound Level: ≤45db

Power Consumption: ≤85w

LCD Display: Present working time, accumulated working time

Alarm: Power Failure alarm, High&low pressure alarm

Net Weight: 6 Kgs

Size: 350×190×320mm



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

Optional

DC12V Car Adapter(1PC)

Lithium Battery (1PC)

Battery charger (1PC)

Trolley bag (1PC)